

Technical Note

Capillary Endfittings



MS Source	Inlet	Outlet
Sciex	1/16" female	1/16" female
Thermo Flex	1/16" female	1/16" female



MS Source	Inlet	Outlet
Bruker CSI	1/16" female	1/16" male
Thermo Easyspray	1/16" female	1/16" male

When purchasing a Fortis capillary the standard configuration is 1/16" female - female fittings. If you require the 1/16" female - male configuration then please request.

Fortis capillaries are available in 75 μ , 150 μ , 300 μ internal diameters. 0.5mm and 1.0mm microbore columns are also available. Capillaries can be packed with any 1.7 μ m, 3 μ m or 5 μ m particle size Fortis, Evosphere or Evosphere BIO stationary phase.

Please see full Capillary Guide for more information.