


Column Evaluation



- Column : **Fortis C18 150x4.6mm 5 μ**
- p/n : **F18-050705**
- Mobile Phase:
60:40 ACN : Water
- Flow : **1ml/min**
- Temp : **25oC**
- λ : **254nm**

1. **Uracil**
2. **Pyrazoline**
3. **Phenol**
4. **Pyridine**
5. **Methylbenzoate**
6. **Dimethylaniline**
7. **Toluene**
8. **Ethylbenzene**
9. **Diethylaniline**